

Hydroxyacyl-coenzyme A dehydrogenase, mitochondrial

Catalog No: #AP84038

Package Size: #AP84038-1 50ug #AP84038-2 100ug #AP84038-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hydroxyacyl-coenzyme A dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HCDH,Medium and short-chain L-3-hydroxyacyl-coenzyme A dehydrogenase,Short-chain 3-hydroxyacyl-CoA dehydrogenase,Had,Hadhsc,Schad
Accession No.	Uniprot ID:Q9WVK7
Uniprot	Q9WVK7
GeneID	113965;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hadh

Note: This product is for in vitro research use only