

Glycogen synthase kinase-3 alpha

Catalog No: #AP84040

Package Size: #AP84040-1 50ug #AP84040-2 100ug #AP84040-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glycogen synthase kinase-3 alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GSK-3 alpha,Factor A,FA,Serine/threonine-protein kinase GSK3A
Accession No.	Uniprot ID:P18265
Uniprot	P18265
GeneID	50686;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gsk3a

Note: This product is for in vitro research use only