

# Insulin-like growth factor-binding protein 1

Catalog No: #AP84041



Package Size: #AP84041-1 50ug #AP84041-2 100ug #AP84041-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Insulin-like growth factor-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IBP-1,Igfbp-1
Accession No.	Uniprot ID:P21743
Uniprot	P21743
GeneID	25685;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Igfbp1

Note: This product is for in vitro research use only