

Interferon regulatory factor 8

Catalog No: #AP84044



Package Size: #AP84044-1 50ug #AP84044-2 100ug #AP84044-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon regulatory factor 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	IRF-8, Interferon consensus sequence-binding protein, ICSBP, ICSBP, ICSBP1
Accession No.	Uniprot ID:Q90871
Uniprot	Q90871
GeneID	396385;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IRF8

Note: This product is for in vitro research use only