

Integrin alpha-IIb

Catalog No: #AP84046



Package Size: #AP84046-1 50ug #AP84046-2 100ug #AP84046-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Integrin alpha-IIb
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GPalpha IIb,GPIIb,Platelet membrane glycoprotein IIb,CD41,GP2B,ITGAB
Accession No.	Uniprot ID:P08514
Uniprot	P08514
GeneID	3674;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITGA2B

Note: This product is for in vitro research use only