

Heme oxygenase 1

Catalog No: #AP84049



Package Size: #AP84049-1 50ug #AP84049-2 100ug #AP84049-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heme oxygenase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HO-1,P32 protein
Accession No.	Uniprot ID:P14901
Uniprot	P14901
GeneID	15368;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmox1

Note: This product is for in vitro research use only