

4-hydroxy-2-oxoglutarate aldolase, mitochondrial

Catalog No: #AP84050



Package Size: #AP84050-1 50ug #AP84050-2 100ug #AP84050-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-hydroxy-2-oxoglutarate aldolase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Dihydrodipicolinate synthase-like,DHDPS-like protein,Probable 2-keto-4-hydroxyglutarate aldolase,Probable KHG-aldolase,Protein 569272,C10orf65,DHDPSL
Accession No.	Uniprot ID:Q86XE5
Uniprot	Q86XE5
GeneID	112817;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HOGA1

Note: This product is for in vitro research use only