

## Interleukin-13 receptor subunit alpha-2

Catalog No: #AP84051



Package Size: #AP84051-1 50ug #AP84051-2 100ug #AP84051-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-13 receptor subunit alpha-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-13 receptor subunit alpha-2, Interleukin-13-binding protein, CD213a2, IL13R
Accession No.	Uniprot ID:Q14627
Uniprot	Q14627
GeneID	3598;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: IL13RA2

Note: This product is for in vitro research use only