

Interleukin-13

Catalog No: #AP84058



Package Size: #AP84058-1 50ug #AP84058-2 100ug #AP84058-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IL-13
Accession No.	Uniprot ID:Q9XSV9
Uniprot	Q9XSV9
GeneID	281247;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL13

Note: This product is for in vitro research use only