

Solute carrier family 2, facilitated glucose transporter member 1

Catalog No: #AP84060

Package Size: #AP84060-1 50ug #AP84060-2 100ug #AP84060-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Solute carrier family 2, facilitated glucose transporter member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Glucose transporter type 1, erythrocyte/brain, GLUT-1, GLUT1
Accession No.	Uniprot ID:P20303
Uniprot	P20303
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC2A1

Note: This product is for in vitro research use only