

Histone deacetylase 1

Catalog No: #AP84062



Package Size: #AP84062-1 50ug #AP84062-2 100ug #AP84062-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone deacetylase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	HD1
Accession No.	Uniprot ID:Q32PJ8
Uniprot	Q32PJ8
GeneID	404126;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HDAC1

Note: This product is for in vitro research use only