

Glutamine synthetase

Catalog No: #AP84066



Package Size: #AP84066-1 50ug #AP84066-2 100ug #AP84066-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glutamine synthetase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	GS, Glutamate decarboxylase, Glutamate--ammonia ligase, p42, GLNS
Accession No.	Uniprot ID:P16580
Uniprot	P16580
GeneID	396489;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GLUL

Note: This product is for in vitro research use only