

Glycogen synthase kinase-3 beta

Catalog No: #AP84068



Package Size: #AP84068-1 50ug #AP84068-2 100ug #AP84068-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glycogen synthase kinase-3 beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GSK-3 beta, Serine/threonine-protein kinase GSK3B
Accession No.	Uniprot ID:Q9WV60
Uniprot	Q9WV60
GeneID	56637;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gsk3b

Note: This product is for in vitro research use only