

## Glutathione S-transferase alpha-3

Catalog No: #AP84069



Package Size: #AP84069-1 50ug #AP84069-2 100ug #AP84069-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glutathione S-transferase alpha-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GST 2-2,GST A3-3,GST AA,Glutathione S-transferase Yc-1,GST Yc1,Gstyc1
Accession No.	Uniprot ID:P04904
Uniprot	P04904
GeneID	108348061;24421;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gsta3

Note: This product is for in vitro research use only