

Indian hedgehog protein

Catalog No: #AP84073



Package Size: #AP84073-1 50ug #AP84073-2 100ug #AP84073-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Indian hedgehog protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IHH,HHG-2
Accession No.	Uniprot ID:P97812
Uniprot	P97812
GeneID	16147;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ihh

Note: This product is for in vitro research use only