

## 4-hydroxy-2-oxoglutarate aldolase, mitochondrial

Catalog No: #AP84080



Package Size: #AP84080-1 50ug #AP84080-2 100ug #AP84080-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	4-hydroxy-2-oxoglutarate aldolase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Dihydrodipicolinate synthase-like,DHDPS-like protein,Probable 2-keto-4-hydroxyglutarate aldolase,Probable KHG-aldolase,DHDPSL
Accession No.	Uniprot ID:Q0P5I5
Uniprot	Q0P5I5
GeneID	506001;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HOGA1

Note: This product is for in vitro research use only