

Glutathione S-transferase P

Catalog No: #AP84084



Package Size: #AP84084-1 50ug #AP84084-2 100ug #AP84084-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutathione S-transferase P
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Chain 7,GST 7-7,GST class-pi
Accession No.	Uniprot ID:P04906
Uniprot	P04906
GeneID	24426;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gstp1

Note: This product is for in vitro research use only