

## Plasma protease C1 inhibitor

Catalog No: #AP84096



Package Size: #AP84096-1 50ug #AP84096-2 100ug #AP84096-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Plasma protease C1 inhibitor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C1 Inh,C1 esterase inhibitor,C1-inhibiting factor,Serpin G1,C1nh
Accession No.	Uniprot ID:P97290
Uniprot	P97290
GeneID	12258;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Serping1

Note: This product is for in vitro research use only