

Intercellular adhesion molecule 4

Catalog No: #AP84102



Package Size: #AP84102-1 50ug #AP84102-2 100ug #AP84102-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Intercellular adhesion molecule 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ICAM-4,CD242
Accession No.	Uniprot ID:Q9ERM2
Uniprot	Q9ERM2
GeneID	78369;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Icam4

Note: This product is for in vitro research use only