

## Beta-hexosaminidase subunit alpha

Catalog No: #AP84109



Package Size: #AP84109-1 50ug #AP84109-2 100ug #AP84109-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Beta-hexosaminidase subunit alpha  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Rt   |
| Other Names        | Beta-N-acetylhexosaminidase subunit alpha, Hexosaminidase subunit A, N-acetyl-beta-glucosaminidase subunit alpha |
| Accession No.      | Uniprot ID:Q641X3  |
| Uniprot            | Q641X3   |
| GeneID             | 300757;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Hexa

Note: This product is for in vitro research use only