

78 kDa glucose-regulated protein

Catalog No: #AP84111



Package Size: #AP84111-1 50ug #AP84111-2 100ug #AP84111-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	78 kDa glucose-regulated protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GRP-78,Heat shock 70 kDa protein 5,Immunoglobulin heavy chain-binding protein,BiP,Steroidogenesis-activator polypeptide,Grp78
Accession No.	Uniprot ID:P06761
Uniprot	P06761
GeneID	25617;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hspa5

Note: This product is for in vitro research use only