

High mobility group protein B1

Catalog No: #AP84115



Package Size: #AP84115-1 50ug #AP84115-2 100ug #AP84115-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | High mobility group protein B1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | High mobility group protein 1,HMG-1 |
| Accession No. | Uniprot ID:Q6YKA4 |
| Uniprot | Q6YKA4 |
| GeneID | 403170; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HMGB1

Note: This product is for in vitro research use only