

## Isopentenyl-diphosphate Delta-isomerase 1

Catalog No: #AP84117



Package Size: #AP84117-1 50ug #AP84117-2 100ug #AP84117-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Isopentenyl-diphosphate Delta-isomerase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Isopentenyl pyrophosphate isomerase 1,IPP isomerase 1
Accession No.	Uniprot ID:P58044
Uniprot	P58044
GeneID	102635781;319554;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Idi1

Note: This product is for in vitro research use only