

Interferon regulatory factor 1

Catalog No: #AP84119

Package Size: #AP84119-1 50ug #AP84119-2 100ug #AP84119-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon regulatory factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IRF-1
Accession No.	Uniprot ID:P10914
Uniprot	P10914
GeneID	3659;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IRF1

Note: This product is for in vitro research use only