

Ubiquitin-like protein ISG15

Catalog No: #AP84121



Package Size: #AP84121-1 50ug #AP84121-2 100ug #AP84121-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin-like protein ISG15
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Interferon-induced 15 kDa protein, Interferon-induced 17 kDa protein, IP17, Ubiquitin cross-reactive protein, hUCRP, G1P2, UCRP
Accession No.	Uniprot ID:P05161
Uniprot	P05161
GeneID	9636;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: ISG15

Note: This product is for in vitro research use only