

## Serine protease inhibitor Kazal-type 7

Catalog No: #AP84123



Package Size: #AP84123-1 50ug #AP84123-2 100ug #AP84123-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine protease inhibitor Kazal-type 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Esophagus cancer-related gene 2 protein, ECRG-2, ECG2, UNQ745/PRO1474
Accession No.	Uniprot ID:P58062
Uniprot	P58062
GeneID	84651;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: SPINK7

Note: This product is for in vitro research use only