

Interferon lambda-3

Catalog No: #AP84125



Package Size: #AP84125-1 50ug #AP84125-2 100ug #AP84125-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon lambda-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IFN-lambda-3, Interleukin-28B, IL-28B, IL28, IL28b
Accession No.	Uniprot ID:Q8CGK6
Uniprot	Q8CGK6
GeneID	338374;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IfnI3

Note: This product is for in vitro research use only