

Interleukin-1 receptor type 2

Catalog No: #AP84126



Package Size: #AP84126-1 50ug #AP84126-2 100ug #AP84126-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-1R-2,CD121 antigen-like family member B,CDw121b,IL-1 type II receptor,Interleukin-1 receptor beta,IL-1R-beta,Interleukin-1 receptor type II,CD121b,IL1RB
Accession No.	Uniprot ID:P27930
Uniprot	P27930
GeneID	7850;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL1R2

Note: This product is for in vitro research use only