

## Interleukin-6 receptor subunit alpha

Catalog No: #AP84127



Package Size: #AP84127-1 50ug #AP84127-2 100ug #AP84127-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                         |
|--------------------|-------------------------------------------------------------------------|
| Product Name       | Interleukin-6 receptor subunit alpha                                    |
| Brief Description  | Recombinant Protein                                                     |
| Host Species       | E.coli                                                                  |
| Purification       | Greater than 90% by SDS-PAGE                                            |
| Species Reactivity | Hu                                                                      |
| Other Names        | IL-6 receptor subunit alpha,IL-6R 1,Membrane glycoprotein 80,gp80,CD126 |
| Accession No.      | Uniprot ID:P08887                                                       |
| Uniprot            | P08887                                                                  |
| GeneID             | 3570;                                                                   |
| Tag Info           | His                                                                     |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                  |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL6R

Note: This product is for in vitro research use only