

Interleukin-1 receptor-like 1

Catalog No: #AP84128



Package Size: #AP84128-1 50ug #AP84128-2 100ug #AP84128-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 receptor-like 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Fit-1,Fos-responsive gene 1 protein,Fit1
Accession No.	Uniprot ID:Q62611
Uniprot	Q62611
GeneID	25556;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Il1rl1

Note: This product is for in vitro research use only