

Inositol monophosphatase 2

Catalog No: #AP84129



Package Size: #AP84129-1 50ug #AP84129-2 100ug #AP84129-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inositol monophosphatase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IMP 2, Inositol-1(or 4)-monophosphatase 2, Myo-inositol monophosphatase A2
Accession No.	Uniprot ID:Q91UZ5
Uniprot	Q91UZ5
GeneID	114663;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Impa2

Note: This product is for in vitro research use only