

## Histone deacetylase 7

Catalog No: #AP84132



Package Size: #AP84132-1 50ug #AP84132-2 100ug #AP84132-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Histone deacetylase 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HD7,Histone deacetylase 7A,HD7a,HDAC7A
Accession No.	Uniprot ID:Q8WUI4
Uniprot	Q8WUI4
GeneID	51564;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HDAC7

Note: This product is for in vitro research use only