

Heat shock 70 kDa protein 1A

Catalog No: #AP84137



Package Size: #AP84137-1 50ug #AP84137-2 100ug #AP84137-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock 70 kDa protein 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Heat shock 70 kDa protein 1,HSP70.1,HSPA1
Accession No.	Uniprot ID:P34930
Uniprot	P34930
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSPA1A

Note: This product is for in vitro research use only