

Interleukin-36 alpha

Catalog No: #AP84145



Package Size: #AP84145-1 50ug #AP84145-2 100ug #AP84145-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-36 alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FIL1 epsilon, Interleukin-1 epsilon, IL-1 epsilon, Interleukin-1 family member 6, IL-1F6, Interleukin-1 homolog 1, IL-1H1, Fil1e, Il1e, Il1f6, Il1h1
Accession No.	Uniprot ID:Q9JLA2
Uniprot	Q9JLA2
GeneID	54448;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL36a

Note: This product is for in vitro research use only