

Heat shock-related 70 kDa protein 2

Catalog No: #AP84146



Package Size: #AP84146-1 50ug #AP84146-2 100ug #AP84146-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock-related 70 kDa protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Heat shock protein 70.2, Testis-specific heat shock protein-related, HST, Hst70
Accession No.	Uniprot ID:P14659
Uniprot	P14659
GeneID	60460;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hspa2

Note: This product is for in vitro research use only