

Interleukin-2 receptor subunit alpha

Catalog No: #AP84148



Package Size: #AP84148-1 50ug #AP84148-2 100ug #AP84148-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-2 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IL-2 receptor subunit alpha,TAC antigen,p55,CD25
Accession No.	Uniprot ID:P12342
Uniprot	P12342
GeneID	281861;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL2RA

Note: This product is for in vitro research use only