

Major histocompatibility complex class I-related gene protein

Catalog No: #AP84151

Package Size: #AP84151-1 50ug #AP84151-2 100ug #AP84151-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Major histocompatibility complex class I-related gene protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MHC class I-related gene protein, Hlals
Accession No.	Uniprot ID:O19477
Uniprot	O19477
GeneID	25119;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mr1

Note: This product is for in vitro research use only