

Interleukin-33

Catalog No: #AP84156



Package Size: #AP84156-1 50ug #AP84156-2 100ug #AP84156-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-33
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	IL-33,Protein DVS27
Accession No.	Uniprot ID:O97863
Uniprot	O97863
GeneID	403810;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL33

Note: This product is for in vitro research use only