

Interleukin-9

Catalog No: #AP84157



Package Size: #AP84157-1 50ug #AP84157-2 100ug #AP84157-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-9,Cytokine P40,T-cell growth factor P40
Accession No.	Uniprot ID:P15247
Uniprot	P15247
GeneID	16198;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL9

Note: This product is for in vitro research use only