

## Interleukin-1 receptor-like 1

Catalog No: #AP84158



Package Size: #AP84158-1 50ug #AP84158-2 100ug #AP84158-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-1 receptor-like 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Interleukin-33 receptor alpha chain,Lymphocyte antigen 84,Protein ST2,Protein T1,Ly84,St2,Ste2
Accession No.	Uniprot ID:P14719
Uniprot	P14719
GeneID	17082;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Il1rl1

Note: This product is for in vitro research use only