

Hemoglobin subunit alpha

Catalog No: #AP84159



Package Size: #AP84159-1 50ug #AP84159-2 100ug #AP84159-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hemoglobin subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Alpha-globin,Hemoglobin alpha chain
Accession No.	Uniprot ID:P01966
Uniprot	P01966
GeneID	100140149;512439;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HBA

Note: This product is for in vitro research use only