

Hydroxyacyl-coenzyme A dehydrogenase, mitochondrial

Catalog No: #AP84161

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84161-1 50ug #AP84161-2 100ug #AP84161-3 1mg

Description

Product Name	Hydroxyacyl-coenzyme A dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	HCDH,Medium and short-chain L-3-hydroxyacyl-coenzyme A dehydrogenase,Short-chain 3-hydroxyacyl-CoA dehydrogenase,HAD,HADHSC,SCHAD
Accession No.	Uniprot ID:P00348
Uniprot	P00348
GeneID	397604;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HADH

Note: This product is for in vitro research use only