

Interleukin-17F

Catalog No: #AP84162



Package Size: #AP84162-1 50ug #AP84162-2 100ug #AP84162-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-17F
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-17F,Cytokine ML-1
Accession No.	Uniprot ID:Q96PD4
Uniprot	Q96PD4
GeneID	112744;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL17F

Note: This product is for in vitro research use only