

Interleukin-6 receptor subunit alpha

Catalog No: #AP84163



Package Size: #AP84163-1 50ug #AP84163-2 100ug #AP84163-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-6 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-6 receptor subunit alpha,IL-6R 1,CD126,Il6r
Accession No.	Uniprot ID:P22272
Uniprot	P22272
GeneID	16194;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Il6ra

Note: This product is for in vitro research use only