

Integrin-linked protein kinase

Catalog No: #AP84164



Package Size: #AP84164-1 50ug #AP84164-2 100ug #AP84164-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin-linked protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	59 kDa serine/threonine-protein kinase,ILK-1,ILK-2,p59ILK,ILK1,ILK2
Accession No.	Uniprot ID:Q13418
Uniprot	Q13418
GeneID	3611;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ILK

Note: This product is for in vitro research use only