

Integrator complex subunit 3

Catalog No: #AP84166



Package Size: #AP84166-1 50ug #AP84166-2 100ug #AP84166-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrator complex subunit 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Int3,SOSS complex subunit A,Sensor of single-strand DNA complex subunit A,SOSS-A
Accession No.	Uniprot ID:Q7TPD0
Uniprot	Q7TPD0
GeneID	229543;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ints3

Note: This product is for in vitro research use only