

Interleukin-1 receptor type 2

Catalog No: #AP84168



Package Size: #AP84168-1 50ug #AP84168-2 100ug #AP84168-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IL-1R-2,CD121 antigen-like family member B,IL-1 type II receptor,Interleukin-1 receptor beta,IL-1R-beta,Interleukin-1 receptor type II,CD121b,IL-1r2,Il1rb
Accession No.	Uniprot ID:P43303
Uniprot	P43303
GeneID	117022;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Il1r2

Note: This product is for in vitro research use only