

Interleukin-1 receptor antagonist protein

Catalog No: #AP84169



Package Size: #AP84169-1 50ug #AP84169-2 100ug #AP84169-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-1 receptor antagonist protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-1RN,IL1 inhibitor,IL-1ra
Accession No.	Uniprot ID:P25085
Uniprot	P25085
GeneID	16181;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Il1rn

Note: This product is for in vitro research use only