

# Interleukin-3 receptor subunit alpha

Catalog No: #AP84170



Package Size: #AP84170-1 50ug #AP84170-2 100ug #AP84170-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-3 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-3 receptor subunit alpha, Interleukin-3 receptor class II alpha chain, CD123, Sut-1
Accession No.	Uniprot ID:P26952
Uniprot	P26952
GeneID	16188;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: IL3ra

Note: This product is for in vitro research use only