

## Interleukin-1 receptor type 2

Catalog No: #AP84172



Package Size: #AP84172-1 50ug #AP84172-2 100ug #AP84172-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-1 receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-1R-2,CD121 antigen-like family member B,IL-1 type II receptor,Interleukin-1 receptor beta,IL-1R-beta,Interleukin-1 receptor type II,CD121b,IL-1r2,Il1rb
Accession No.	Uniprot ID:P27931
Uniprot	P27931
GeneID	16178;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Il1r2

Note: This product is for in vitro research use only